



IVTScrip™ mRNA-Human AKR1D1, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18378MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18378MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the AKR1D1 protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001190906.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001190906.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CBAS2; SRD5B1; 3o5bred
Description	The enzyme encoded by this gene is responsible for the catalysis of the 5-beta-reduction of bile acid intermediates and steroid hormones carrying a delta(4)-3-one structure. Deficiency of this enzyme may contribute to hepatic dysfunction. Three transcript variants encoding different isoforms have been found for this gene. Other variants may be present, but their full-length natures have not been determined yet. [provided by RefSeq, Jul 2010]