



IVTScrip™ mRNA-Human KRT19, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK1840MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1840MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the KRT19 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002276
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002276
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CK19; MGC1536
Description	This gene is a protein coding gene. A structural constituent of muscle. Involved in response to estr ogen and sarcomere organization. Localizes to costamere and plasma membrane.