



IVTScrip™ mRNA-Human AMER3, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18449MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18449MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the AMER3 protein. This product can be used in Secretory progenitor cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001105193.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001105193.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	FAM123C
Description	AMER3 (APC Membrane Recruitment Protein 3) is a Protein Coding gene. Diseases associated with AMER3 include Tracheoesophageal Fistula With Or Without Esophageal Atresia and Osteopathia Striata With Cranial Sclerosis. Gene Ontology (GO) annotations related to this gene include lipid binding. An important paralog of this gene is AMER1.