




IVTScrip™ mRNA-Human CD11a, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK1844MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1844MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the CD11a protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002209
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002209
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CD11A; LFA1A
Description	This gene is a protein coding gene. Exhibits ICAM-3 receptor activity and cell adhesion molecule binding activity. Involved in several processes, including T cell activation via T cell receptor contact with antigen bound to MHC molecule on antigen presenting cell; heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules; and memory T cell extravasation. Localizes to cell surface and integrin alphaL-beta2 complex.