



IVTScrip™ mRNA-Human AMY2A, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18462MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18462MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It encodes the AMY2A protein. This product can be used in Acinar cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000699.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_000699.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	PA; AMY2
Description	This gene encodes a member of the alpha-amylase family of proteins. Amylases are secreted proteins that hydrolyze 1,4-alpha-glucoside bonds in oligosaccharides and polysaccharides, catalyzing the first step in digestion of dietary starch and glycogen. This gene and several family members are present in a gene cluster on chromosome 1. This gene encodes an amylase isoenzyme produced by the pancreas. [provided by RefSeq, Jan 2015]