




IVTScrip™ mRNA-Human CD195, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK1848MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1848MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the CD195 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_0017705
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_0017705
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CMKBR5; CD195
Description	This gene is a protein coding gene. Exhibits several functions, including C-C chemokine receptor activity; chemokine (C-C motif) ligand 5 binding activity; and identical protein binding activity. Involved in several processes, including cellular response to lipopolysaccharide; intracellular signal transduction; and release of sequestered calcium ion into cytosol by sarcoplasmic reticulum. Localizes to endosome and external side of plasma membrane. Implicated in several diseases, including Kawasaki disease; West Nile fever; autoimmune disease (multiple); berylliosis; and sarcoidosis (multiple). Biomarker of several diseases, including autoimmune disease (multiple); dry eye syndrome (multiple); lung disease (multiple); pleural tuberculosis; and respiratory syncytial virus infectious disease.