




IVTScrip™ mRNA-Human ANKUB1, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18841MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18841MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the ANKUB1 protein. This product can be used in Ciliated epithelial cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001144960.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001144960.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	C3orf16
Description	ANKUB1 (Ankyrin Repeat And Ubiquitin Domain Containing 1) is a Protein Coding gene. An important paralog of this gene is ZFAND4.