



IVTScrip™ mRNA-Human APBB1IP, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK19199MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK19199MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the APBB1IP protein. This product can be used in Microglial cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_019043.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_019043.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	RIAM; INAG1; PREL1; RARP1
Description	APBB1IP (Amyloid Beta Precursor Protein Binding Family B Member 1 Interacting Protein) is a Protein Coding gene. Diseases associated with APBB1IP include Leukocyte Adhesion Deficiency, Type I and Cerebral Cavernous Malformations. Among its related pathways are Signaling by GPCR and Response to elevated platelet cytosolic Ca ²⁺ . An important paralog of this gene is RAPH1.