



IVTScrip™ mRNA-Human ANKRD26P1, (Cap 1, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK19411MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK19411MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the ANKRD26P1 protein. This product can be used in Gonadal mitotic phase fetal germ cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	FLJ43980
Description	ANKRD26P1 (Ankyrin Repeat Domain 26 Pseudogene 1) is a Pseudogene.