



IVTScrip™ mRNA-Human AKIRIN2, (Cap 1, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK19570MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK19570MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP & Pseudo-UTP. It encodes the AKIRIN2 protein. This product can be used in Regulatory T (Treg) cell-related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_018064.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_018064.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	FBI1; C6orf166; dJ486L4.2
Description	AKIRIN2 (Akirin 2) is a Protein Coding gene. Diseases associated with AKIRIN2 include Tick Infestation and Metacarpal 4-5 Fusion. Gene Ontology (GO) annotations related to this gene include enzyme binding. An important paralog of this gene is AKIRIN1.