



IVTScrip™ mRNA-Human AMACR, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK1974MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1974MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the AMACR protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001167595.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001167595.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	CBAS4; P504S; 2-methylacyl-CoA racemase; AMACRD; RACE; RM
Description	This gene is a protein coding gene. Exhibits alpha-methylacyl-CoA racemase activity and signaling receptor binding activity. Involved in bile acid metabolic process. Localizes to mitochondrion; peroxisome; and plasma membrane. Implicated in alpha-methylacyl-CoA racemase deficiency; congenital bile acid synthesis defect 4; and urinary bladder cancer. Biomarker of breast carcinoma; prostate carcinoma in situ; and renal cell carcinoma.