



IVTScrip™ mRNA-Human AIM1, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK19833MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK19833MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the AIM1 protein. This product can be used in Oogenesis phase fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001371242.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001371242.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ST4; CRYBG1
Description	CRYBG1 (Crystallin Beta-Gamma Domain Containing 1) is a Protein Coding gene. Diseases associated with CRYBG1 include Melanoma. An important paralog of this gene is CRYBG3.