



IVTScrip™ mRNA-Human ANKMY1, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK19994MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK19994MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ANKMY1 protein. This product can be used in Ciliated epithelial cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001282771.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001282771.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ZMYND13
Description	ANKMY1 (Ankyrin Repeat And MYND Domain Containing 1) is a Protein Coding gene. Diseases associated with ANKMY1 include Venous Hemangioma and Sleeping Sickness.