



IVTScrip™ mRNA-Human ANKRD49, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20029MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20029MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ANKRD49 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_017704.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_017704.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	FGIF; GBIF
Description	ANKRD49 (Ankyrin Repeat Domain 49) is a Protein Coding gene. Diseases associated with ANKR D49 include Myasthenic Syndrome, Congenital, 10 and Syndromic X-Linked Intellectual Disability Snyder Type.