



IVTScrip™ mRNA-Human ANKRD6, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20032MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20032MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ANKRD6 protein. This product can be used in Trophectoderm cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001242809.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001242809.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKRD6 (Ankyrin Repeat Domain 6) is a Protein Coding gene. Diseases associated with ANKRD6 include Neural Tube Defects and Ocular Motility Disease. An important paralog of this gene is ANKRD10.
--------------------	---