




IVTScrip™ mRNA-Human ANKRD66, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20035MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20035MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ANKRD66 protein. This product can be used in Ciliated cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001162435.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001162435.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKRD66 (Ankyrin Repeat Domain 66) is a Protein Coding gene.
--------------------	--