



IVTScrip™ mRNA-Human AJUBA, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20140MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20140MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It codes the AJUBA protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001289097.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001289097.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	JUB
Description	AJUBA (Ajuba LIM Protein) is a Protein Coding gene. Diseases associated with AJUBA include Wil son-Turner X-Linked Mental Retardation Syndrome and Bjornstad Syndrome. Among its related pathways are Hippo signaling pathway and Regulation of PLK1 Activity at G2/M Transition. Gene Ontology (GO) annotations related to this gene include chromatin binding and actin filament binding. A n important paralog of this gene is WTIP.