



IVTScrip™ mRNA-Human ALG8, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20218MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20218MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & N1-Methylpseudo-UTP. It codes the ALG8 protein. This product can be used in Liver bud hepatic cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001007027.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001007027.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CDG1H; PCLD3
Description	This gene encodes a member of the ALG6/ALG8 glucosyltransferase family. The encoded protein catalyzes the addition of the second glucose residue to the lipid-linked oligosaccharide precursor for N-linked glycosylation of proteins. Mutations in this gene have been associated with congenital disorder of glycosylation type 1h (CDG-1h). Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]