



IVTScrip™ mRNA-Human HGAL, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2034MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2034MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the HGAL protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001190259.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001190259.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Description	This gene is a protein coding gene. Exhibits actin binding activity; myosin II binding activity; and protein kinase binding activity. Involved in negative regulation of lymphocyte migration and regulation of B cell receptor signaling pathway. Predicted to localize to cytoplasm and plasma membrane.
--------------------	--