



IVTScrip™ mRNA-Human AGMO, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20404MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20404MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the AGMO protein. This product can be used in Mitotic arrest phase fetal germ cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001004320.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001004320.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	TMEM195
Description	The protein encoded by this gene is a tetrahydrobiopterin- and iron-dependent enzyme that cleaves the ether bond of alkylglycerols. Sequence comparisons distinguish this protein as forming a third, distinct class of tetrahydrobiopterin-dependent enzymes. Variations in this gene have been associated with decreased glucose-stimulated insulin response, type 2 diabetes, and susceptibility to intracranial aneurysms. [provided by RefSeq, Aug 2012]