



IVTScrip™ mRNA-Human AIMP2, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20436MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20436MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the AIMP2 protein. This product can be used in Liver and hepatic cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001326606.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001326606.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	P38; JTV1; HLD17; JTV-1
Description	The protein encoded by this gene is part of the aminoacyl-tRNA synthetase complex, which contains nine different aminoacyl-tRNA synthetases and three non-enzymatic factors. The encoded protein is one of the non-enzymatic factors and is required for assembly and stability of the complex. [provided by RefSeq, May 2016]