



IVTScrip™ mRNA-Human ANGEL1, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20577MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20577MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methoxy-UTP. It encodes the ANGEL1 protein. This product can be used in Leydig cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001370746.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001370746.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	Ccr4e; KIAA0759
Description	ANGEL1 (Angel Homolog 1) is a Protein Coding gene. Diseases associated with ANGEL1 include Arrhythmogenic Right Ventricular Dysplasia, Familial, 8 and Juvenile Glaucoma. An important para log of this gene is ANGEL2.