




IVTScrip™ mRNA-Human AHSA1, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK20726MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK20726MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 2-Thio-UTP. It encodes the AHSA1 protein. This product can be used in Ciliated cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001321441.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001321441.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	p38; AHA1; hAha1; C14orf3
Description	AHSA1 (Activator Of HSP90 ATPase Activity 1) is a Protein Coding gene. Diseases associated with AHSA1 include Stomatitis and Cystic Fibrosis. Among its related pathways are Mechanisms of C FTR activation by S-nitrosoglutathione (normal and CF) and T cell receptor signaling pathway. Gene Ontology (GO) annotations related to this gene include chaperone binding and ATPase activator activity.