




IVTScrip™ mRNA-Human OGN, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2091MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2091MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the OGN protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_033014
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_033014
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	SLRR3A
Description	This gene is a protein coding gene. Predicted to have extracellular matrix structural constituent conferring compression resistance and growth factor activity. Predicted to be involved in negative regulation of smooth muscle cell proliferation. Localizes to collagen-containing extracellular matrix and extracellular exosome.