




IVTScrip™ mRNA-Human AGTPBP1, (Cap 1, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK21017MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK21017MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the AGTPBP1 protein. This product can be used in Effector CD8+ memory T (Tem) cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001286715.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001286715.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CCP1; NNA1; CONDCA
Description	NNA1 is a zinc carboxypeptidase that contains nuclear localization signals and an ATP/GTP-binding motif that was initially cloned from regenerating spinal cord neurons of the mouse.[supplied by O MIM, Jul 2002]