



IVTScrip™ mRNA-Human ANP32E, (Cap 1, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK21252MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK21252MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the ANP32E protein. This product can be used in Neural progenitor cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001136478.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001136478.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	LANPL; LANP-L
Description	ANP32E (Acidic Nuclear Phosphoprotein 32 Family Member E) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include histone binding and phosphatase inhibitor activity. An important paralog of this gene is ANP32A.