



## IVTScrip™ mRNA-Human AP1AR, (Cap 1, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK21272MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK21272MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 1 & 5-Methyl-CTP. It encodes the AP1AR protein. This product can be used in Mitotic arrest phase fetal germ cell-related researches.

### Specifications

<b>Modified bases</b>	5-Methyl-CTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001128426.3
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

### SPECIFICATIONS

<b>Modified bases</b>	5-Methyl-CTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001128426.3
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

## GENE INFORMATION

<b>Alternative Names</b>	2C18; GBAR; C4orf16; PRO0971; gamma-BAR
<b>Description</b>	AP1AR (Adaptor Related Protein Complex 1 Associated Regulatory Protein) is a Protein Coding gene. Diseases associated with AP1AR include Gliosarcoma. Gene Ontology (GO) annotations related to this gene include kinesin binding and AP-1 adaptor complex binding.