



## IVTScrip™ mRNA-Human APOO, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK21639MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK21639MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the APOO protein. This product can be used in Ciliated cell-related researches.

#### Specifications

|                |                       |
|----------------|-----------------------|
| Modified bases | Pseudo-UTP            |
| 5' Cap         | Cap 0                 |
| Species        | Human                 |
| RefSeq         | NM_024122.5           |
| Applications   | Gene therapy research |
| Format         | Powder                |
| Quantity       | 100 µg                |
| Purification   | Chromatography        |

#### SPECIFICATIONS

|                       |                       |
|-----------------------|-----------------------|
| <b>Modified bases</b> | Pseudo-UTP            |
| <b>5' Cap</b>         | Cap 0                 |
| <b>Species</b>        | Human                 |
| <b>RefSeq</b>         | NM_024122.5           |
| <b>Applications</b>   | Gene therapy research |
| <b>Format</b>         | Powder                |
| <b>Quantity</b>       | 100 µg                |
| <b>Purification</b>   | Chromatography        |

## GENE INFORMATION

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | MIC26; Mic23; My025; FAM121B; MICOS26                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Description</b>       | This gene is a member of the apolipoprotein family. Members of this protein family are involved in the transport and metabolism of lipids. The encoded protein associates with HDL, LDL and VLDL lipoproteins and is characterized by chondroitin-sulfate glycosylation. This protein may be involved in preventing lipid accumulation in the myocardium in obese and diabetic patients. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 3, 4, 5, 12 and 16.[provided by RefSeq, Sep 2009] |