



IVTScrip™ mRNA-Human ARC, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK21659MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK21659MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ARC protein. This product can be used in Oligodendrocyte progenitor cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_015193.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_015193.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description ARC (Activity Regulated Cytoskeleton Associated Protein) is a Protein Coding gene. Diseases associated with ARC include Hepatocellular Carcinoma and Gordon Holmes Syndrome. Among its related pathways are Amphetamine addiction and BDNF-TrkB Signaling.