




IVTScrip™ mRNA-Human CDW198, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2166MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2166MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the CDW198 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_005201
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_005201
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	GPRCY6; CKR-L1; CKRL1; CY6
Description	This gene is a protein coding gene. Predicted to have C-C chemokine binding activity and C-C chemokine receptor activity. Predicted to be involved in several processes, including calcium-mediated signaling; cell chemotaxis; and positive regulation of cytosolic calcium ion concentration. Predicted to localize to external side of plasma membrane.