



IVTScrip™ mRNA-Human CD13, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2194MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2194MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the CD13 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001150
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001150
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	PEPN; CD13
Description	This gene is a protein coding gene. Predicted to have metalloaminopeptidase activity; peptide binding activity; and zinc ion binding activity. Predicted to be involved in several processes, including peptide catabolic process; proteolysis; and regulation of blood pressure. Localizes to endoplasmic reticulum-Golgi intermediate compartment; extracellular exosome; and plasma membrane. Biomarker of severe acute respiratory syndrome.