



IVTScrip™ mRNA-Human ARFGAP2, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK21970MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK21970MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ARFGAP2 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001242832.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001242832.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	IRZ; ZFP289; ZNF289; NBLA10535
Description	ARFGAP2 (ADP Ribosylation Factor GTPase Activating Protein 2) is a Protein Coding gene. Diseases associated with ARFGAP2 include Hypertrophic Cardiomyopathy and Autoimmune Lymphoproliferative Syndrome. Among its related pathways are Vesicle-mediated transport and Metabolism of proteins. Gene Ontology (GO) annotations related to this gene include GTPase activator activity. An important paralog of this gene is ARFGAP3.