



IVTScrip™ mRNA-Human ARHGAP45, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK22000MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK22000MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ARHGAP45 protein. This product can be used in Naive CD8+ T cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001258328.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001258328.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGAP45 (Rho GTPase Activating Protein 45) is a Protein Coding gene. Diseases associated with ARHGAP45 include Graft-Versus-Host Disease and Acute Graft Versus Host Disease. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling. An important paralog of this gene is ARHGAP29.
--------------------	--