




IVTScrip™ mRNA-Human ARL6IP6, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK22050MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK22050MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ARL6IP6 protein. This product can be used in Plasmacytoid dendritic cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001350068.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001350068.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	AIP6; AIP-6; PFAAP1
Description	ARL6IP6 (ADP Ribosylation Factor Like GTPase 6 Interacting Protein 6) is a Protein Coding gene. Diseases associated with ARL6IP6 include Cutis Marmorata Telangiectatica Congenita and Coffin-Siris Syndrome 2.