



IVTScrip™ mRNA-Human ASZ1, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK22134MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK22134MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ASZ1 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001301821.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001301821.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ALP1; GASZ; Orf3; ANKL1; C7orf7; CT1.19
Description	ASZ1 (Ankyrin Repeat, SAM And Basic Leucine Zipper Domain Containing 1) is a Protein Coding gene. Diseases associated with ASZ1 include Hepatitis E and Leukodystrophy, Hypomyelinating, 1. 2. Among its related pathways are Mitotic Prophase and PIWI-interacting RNA (piRNA) biogenesis. Gene Ontology (GO) annotations related to this gene include obsolete signal transducer activity. An important paralog of this gene is ANKRD34C.