




IVTScrip™ mRNA-Human ATAD5, (Cap 0, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK22141MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK22141MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ATAD5 protein. This product can be used in Goblet cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_024857.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_024857.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ELG1; FRAG1; C17orf41
Description	ATAD5 (ATPase Family AAA Domain Containing 5) is a Protein Coding gene. Diseases associated with ATAD5 include Chromosome 17Q11.2 Deletion Syndrome, 1.4-Mb and Cafe-Au-Lait Spots, Multiple. An important paralog of this gene is CRLF3.