




IVTScrip™ mRNA-Human IL7R, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2249MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2249MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the IL7R protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002185
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002185
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	IL-7RA; CD127; CDW127; IL7RA
Description	This gene is a protein coding gene. Exhibits cytokine receptor activity. Involved in positive regulation of cell population proliferation and positive regulation of receptor signaling pathway via STAT. Predicted to localize to clathrin-coated vesicle membrane and plasma membrane. Implicated in severe combined immunodeficiency and severe combined immunodeficiency, autosomal recessive, T cell-negative, B cell-positive, Nk cell-positive.