



IVTScrip™ mRNA-Human ARF5, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK22867MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK22867MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ARF5 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001662.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001662.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	This gene is a member of the human ADP-ribosylation factor (ARF) gene family. These genes encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The gene products include 6 ARF proteins and 11 ARF-like proteins and constitute 1 family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2, and ARF3), class II (ARF4 and ARF5) and class III (ARF6). The members of each class share a common gene organization. [provided by RefSeq, Dec 2010]
--------------------	--