



IVTScrip™ mRNA-Human ARHGEF33, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK22917MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK22917MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ARHGEF33 protein. This product can be used in Oogenesis phase fetal germ cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001145451.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001145451.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGEF33 (Rho Guanine Nucleotide Exchange Factor 33) is a Protein Coding gene. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling. Gene Ontology (GO) annotations related to this gene include Rho guanyl-nucleotide exchange factor activity.
--------------------	---