




IVTScrip™ mRNA-Human ASUN, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK23030MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK23030MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ASUN protein. This product can be used in Liver and hepatic cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_018164.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_018164.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	INTS13; GCT1; NET48; Mat89Bb; SPATA30; C12orf11
Description	INTS13 (Integrator Complex Subunit 13) is a Protein Coding gene. Diseases associated with INTS13 include Juvenile Type Testicular Granulosa Cell Tumor and Testicular Granulosa Cell Tumor. Among its related pathways are Formation of HIV-1 elongation complex containing HIV-1 Tat and Gene Expression.