



IVTScrip™ mRNA-Human ATAD3C, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK23040MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK23040MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ATAD3C protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001039211.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001039211.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ATAD3C (ATPase Family AAA Domain Containing 3C) is a Protein Coding gene. Diseases associated with ATAD3C include Pontocerebellar Hypoplasia, Type 4 and Pontocerebellar Hypoplasia, Type 3. Gene Ontology (GO) annotations related to this gene include four-way junction helicase activity. An important paralog of this gene is ATAD3A.
--------------------	--