




IVTScrip™ mRNA-Human RPE65, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2314MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2314MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the RPE65 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000329
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000329
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	mRPE65; BCO3; RP20
Description	This gene is a protein coding gene. Exhibits all-trans-retinyl-palmitate hydrolase, 11-cis retinol forming activity and isomerase activity. Involved in detection of light stimulus involved in vision; retina homeostasis; and terpenoid metabolic process. Localizes to endoplasmic reticulum membrane. Implicated in Leber congenital amaurosis 2; retinitis pigmentosa; and retinitis pigmentosa 20. Biomarker of basal cell carcinoma and squamous cell carcinoma.