



## IVTScrip™ mRNA-Human ARHGAP20, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK23186MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK23186MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ARHGAP20 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

### Specifications

|                |                           |
|----------------|---------------------------|
| Modified bases | 5-Methyl-CTP & Pseudo-UTP |
| 5' Cap         | Cap 0                     |
| Species        | Human                     |
| RefSeq         | NM_001258415.2            |
| Applications   | Gene therapy research     |
| Format         | Powder                    |
| Quantity       | 100 µg                    |
| Purification   | Chromatography            |

### SPECIFICATIONS

|                       |                           |
|-----------------------|---------------------------|
| <b>Modified bases</b> | 5-Methyl-CTP & Pseudo-UTP |
| <b>5' Cap</b>         | Cap 0                     |
| <b>Species</b>        | Human                     |
| <b>RefSeq</b>         | NM_001258415.2            |
| <b>Applications</b>   | Gene therapy research     |
| <b>Format</b>         | Powder                    |
| <b>Quantity</b>       | 100 µg                    |
| <b>Purification</b>   | Chromatography            |

## GENE INFORMATION

|                          |                                                                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | RARHOGAP                                                                                                                                                                    |
| <b>Description</b>       | The protein encoded by this gene is an activator of RHO-type GTPases, transducing a signal from RAP1 to RHO and impacting neurite outgrowth. [provided by RefSeq, Sep 2016] |