



IVTScrip™ mRNA-Human CD27, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2319MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2319MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the CD27 protein. This product can be used in -related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001242
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001242
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	T14; MGC20393; Tp55; TNFRSF7
Description	This gene is a protein coding gene. Exhibits cysteine-type endopeptidase inhibitor activity involved in apoptotic process. Involved in extrinsic apoptotic signaling pathway and negative regulation of a apoptotic process. Predicted to localize to external side of plasma membrane. Implicated in lympho proliferative syndrome 2.