



IVTScrip™ mRNA-Human ATG2A, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK23657MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK23657MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ATG2A protein. This product can be used in Mitotic fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001367971.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001367971.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ATG2A (Autophagy Related 2A) is a Protein Coding gene. Diseases associated with ATG2A include Auditory Agnosia. Among its related pathways are Pathways of neurodegeneration - multiple diseases and Autophagy Pathway. An important paralog of this gene is ATG2B.
--------------------	---