



IVTScrip™ mRNA-Human ARHGAP39, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK24398MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK24398MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ARHGAP39 protein. This product can be used in Ciliated cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001308207.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001308207.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ARHGAP39 (Rho GTPase Activating Protein 39) is a Protein Coding gene. Diseases associated with ARHGAP39 include Charcot-Marie-Tooth Disease, X-Linked Recessive, 3 and Neuronopathy, Distal Hereditary Motor, Type I. Among its related pathways are Signaling by GPCR and p75 NTR receptor-mediated signalling.
--------------------	--