




IVTScrip™ mRNA-Human ARL5A, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK24444MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK24444MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ARL5A protein. This product can be used in CD4+ cytotoxic T cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001037174.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001037174.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ARL5; ARFLP5
Description	The protein encoded by this gene belongs to the ARF family of GTP-binding proteins. With its distinctive nuclear/nucleolar localization and interaction with HP1alpha, the protein is developmentally regulated and may play a role(s) in nuclear dynamics and/or signaling cascades during embryonic development. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene has multiple pseudogenes. [provided by RefSeq, Jul 2008]