



IVTScrip™ mRNA-Human CD209, (Cap 0, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2481MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK2481MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the CD209 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_1528669
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_1528669
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	DC-SIGN; CLEC4L
Description	This gene is a protein coding gene. Predicted to have mannose binding activity and virion binding activity. Involved in B cell adhesion and positive regulation of T cell proliferation. Localizes to external side of plasma membrane. Implicated in Coronavirus infectious disease; aspergillosis; dengue disease; human immunodeficiency virus infectious disease; and tuberculosis. Biomarker of asthma; pulmonary hypertension; and tuberculosis.