



## IVTScrip™ mRNA-Human CD66a, (Cap 0, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK2514MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK2514MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the CD66a protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	5-Methoxy-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001024912
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

#### SPECIFICATIONS

<b>Modified bases</b>	5-Methoxy-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001024912
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

<b>Alternative Names</b>	CEACAM1; BGP; BGP1; BGPI
<b>Description</b>	This gene is a protein coding gene. This gene encodes a member of the carcinoembryonic antigen (CEA) gene family, which belongs to the immunoglobulin superfamily. Two subgroups of the CEA family, the CEA cell adhesion molecules and the pregnancy-specific glycoproteins, are located within a 1.2 Mb cluster on the long arm of chromosome 19.