




IVTScrip™ mRNA-Human ATP1A1-AS1, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK25183MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK25183MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ATP1A1-AS1 protein. This product can be used in FGFR1HighNME5- epithelial cell-related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ATP1A1OS; C1orf203
Description	ATP1A1-AS1 (ATP1A1 Antisense RNA 1) is an RNA Gene, and is affiliated with the lncRNA class. Diseases associated with ATP1A1-AS1 include Charcot-Marie-Tooth Disease, Axonal, Type 2Dd and Hypomagnesemia, Seizures, And Mental Retardation 2.